
Reactiv / Kit Masson thricrome 100 teste, produs original Bio-Optica/Italia, certificare IVD

I.D.: 94710375

Data publicarii 14.03.24

Coduri CPV

24224000-0

Pretul estimativ: 487,00 RON - 487,00 RON

Descriere: Kit Masson tricromica cu Anilina blue(04-010802)sau verde deschis(04-011802) Kit pentru 100 de teste,timp de realizare 35 de minute; Garantie 2 ani. Intracellular fibrils, negative image of Golgi apparatus. It is very suitable for black and white microphotographs. Bring section to distilled water. Put on the section 6 drops of reagent A and add 6 drops of reagent B: leave to act 10 minutes. Without washing, drain the slide and put on the section 10 drops of reagent C: leave to act 4 minutes. Wash quickly (3-4 seconds) in distilled water and put on the section 10 drops of reagent D, leave to act 4 minutes. Wash in distilled water and put on the section 10 drops of reagent E: leave to act 10 minutes. Without washing, drain the slide and put 10 drops of reagent F: leave to act 5 minutes. Wash in distilled water and dehydrate rapidly through ascending alcohols, stop for 1 minute at the last absolute ethanol. Clear in xylene and mount. Rezultate Nuclei și gameți.....negru Citoplasmă, keratină, fibre musculare, granule acidofile....rosu Colagen, mucus, granule bazofile ale hipofizei...verde Granule de celule Delta hipofizei.....verde eritrocite.....galben Reactivi: Soluția de fier hematoxină Weigert - B 18 ml Hematoxilina de fier Weigert - O soluție 18 ml Solutie stabila alcoolica cu acid picric 30 ml Acidul Ponceau fuchsin conform Masson 30 ml Soluție de acid fosfomolbdic 30 ml Soluție verde deschis conform Goldner 30 ml References Bio-Optica/Italia Podusul (Dispozitivul medical)
